





Directions:

Do not dilute. Dampen a cleanroomapproved, lint-free cloth, sponge, or mop with the 5002 cleaner and wipe hard surfaces or articles. Allow to dry thoroughly. Do not soak or wet the surfaces; use only enough to cause a slight moisture layer.

Note:

This product is not cleanroom packaged in double bags.

SDS available at www.aclstaticide.com.

Product# 5002: 1-quart bottle; 12 per case



5002 Staticide Clean Room®

Anti-static cleaning for critical environments

Staticide® Clean Room is a specialized anti-static cleaner designed for cleanliness in sensitive environments. It removes static electricity from floors, chairs, work surfaces, glass, and plastics while reducing static caused by friction. By preventing static attraction to dust and dirt, ACL's purest static control topical helps to keep surfaces cleaner longer.

Ideal for the following applications:

- Electronic manufacturing facilities
- Laboratories
- Healthcare

FEATURES

- **Filtered to less than 0.5 microns**
- Product is not intended for food contact surfaces or equipment
- **Non-corrosive formula is chloride- and sodium-free**
- **Suitable for industrial Class 100 (ISO Class 5) cleanrooms**
- Low-charging, anti-static formula
- ***** Color is translucent amber
- **B** Decays static in less than 2 seconds



